

### Sterilizer Monitoring Log - Quality Control

**Facility** \_\_\_\_\_  
**Sterilizer ID** \_\_\_\_\_  
**Model & Serial #** \_\_\_\_\_

Date	Initials	Load #	Items	Temp °C	Pressure	Cycle Length	Cycle Completed		Test Strips Used		Test Results Achieved		Comments/ Corrective Actions
						Min.	Yes✓	No✓	Yes✓	No✓	Yes✓	No✓	

**Note:** If the sterilizer monitoring log is kept beside the sterilizer and completed for each load, it should serve as your record of sterilization. Please check with your Public Health Unit regarding your reporting requirements.

### **Explanation of Headings**

**Date:** Record the date of sterilization

**Initials:** Record the initials of the person loading and running the sterilizer

**Load #:** Record the # of the load

**Items:** List the items being sterilized

**Temperature °C:** Record the temperature in °C

**Pressure:** Record the pressure

**Cycle Length:** Record the cycle length in minutes

**Cycle Completed:** Check yes if the cycle was completed without interruption

**Test Strips Used:** Check yes if the test strips were placed in the sterilizer

**Test Results Achieved:** Once you have received the results from Near North Laboratories, check yes if the strips passed

**Comments/Corrective Actions:**